

Spirometry simplified

We are your primary partner in spirometry. We know how vital your work is and that real people with real needs are reliant on your care. It's why we strive to make your work and your patients' health journeys easier.



Thank you for choosing us!

Call: +358 44 7777 002
[medikro.com](https://www.medikro.com)



Medikro Oy
Pioneerinkatu 3
FI-70800 Kuopio
Finland

ISO
13485

M8178-1.0-en-gb, Medikro Accessories, product leaflet
2024-09-18

Spirometry Accessories

Complete your spirometry system

MEDIKRO | Spirometry
simplified

Our accessories



SpiroSafe

Product code: M9256
Disposable flow transducer



Pressure Tube

Product code: M9264-200
Stabilize the test signal



Grip

Product code: M91227-25
Makes every measurement easier



Calibration Syringe

Product code: M9474
High-precision volume calibration syringe: volume 3000ml



Bacterial & Viral Filter

Product code: M9293
Protection for Healthcare personnel



Nose Clip

Product code: MRN-03
Prevent test-quality loss



Ordering packages

Streamline your workflow with these ready-to-use study kits.



Spirometry Study Kit 01

- 100pcs of M9256 SpiroSafe Flow Transducers
- 100pcs of M9293 Bacterial & Viral Filters
- 100pcs of MRN Nose Clips



Spirometry Study Kit 02

- 100pcs of M9256 SpiroSafe Flow Transducers
- 100pcs of M9293 Bacterial & Viral Filters

Our spirometers

Screening

- Medikro Duo is a 2-in-1 device for asthma and COPD screening.

Diagnostic

- Medikro Pro is an accurate and precise diagnostic spirometer with embedded ambient condition measuring sensors for automatic real-time BTPS correction
- Medikro Primo is a professional and efficient PC-based spirometer for primary care.

