

Medikro Duo

PRODUCT CODE
M920

The 2-in-1 portable spirometer for asthma
and COPD screening

Quick test mode

Test quality indicators

Calibration free

User-Friendly





PURCHASE DETAILS

Product code – M920


Medikro Duo Spirometer with ambient conditions measurement sensors

- M9356-180W USB Cable C Plug to A Plug, length 1.8m, color white
- M9831 Medikro Spirometry Software
- M8304 User's Manual
- M8141-1.0 Installation Guide

Medikro Duo is a versatile device that can be used for both asthma and COPD screening. With its quick test mode, spirometry screening tests can be completed in just two minutes, and Medikro's quality indicators ensure both speed and accuracy.

Medikro Duo is utilizing Medikro Spirometry Software (MSS). Additionally, MSS makes possible to view and monitor screening data from multiple spirometers in one single interface. The Medikro Duo requires no calibration, and utilises the reliable and disposable Medikro SpiroSafe flow transducers.

HIGHLIGHTS



1 3 measurement modes

The Medikro Duo is a versatile spirometer that offers three different operation modes to choose from: FVC, FEV6 and PEF. This flexibility makes it suitable for a wide range of screening purposes.

2 Quick test mode

The Medikro Duo is a handy tool that lets you perform spirometry screening in less than two minutes, thanks to its separate quick test mode. This feature allows you to complete screening quickly and efficiently, making it a great option for busy clinics or non-clinical settings.

3 Screening quality indicators

These indicators can detect hesitation, lazy blowing, coughing, and early termination during spirometry screening, providing you with valuable insights into the quality of the results.

SERVICE AGREEMENT

Opt for a service agreement with Medikro to enhance your equipment's longevity and reliability.

Indefinite Warranty

Service agreement customers receive a warranty extension as long as they continue annual maintenance, ensuring long-term value.

Hassle-Free Maintenance

We handle scheduling and reminders, allowing you to focus on core responsibilities while maintaining optimal device performance.

CUSTOMERS' FEEDBACK

“Medikro systems are reliable and affordable, their measurements precise and the software well thought out!”

ENRIQUE CIMAS. GP FAMILY MEDICINE



Medikro Oy
Pioneerinkatu 3
70800 Kuopio, Finland

+358 44 777 7002
www.medikro.com